

Electrical Characterization of Semiconductor Nanostructures for Spintronics Applications



Intern : Jan Rehorik

Major : Computer Engineering

Mentor : Jason Stephens

Faculty Advisor : David Awschalom



Funded By : Defense Advanced Research Projects Agency (DARPA)

Semiconductors



- SC resistivity between that of conductors and insulators
- Resistivity can be tailored over many orders of magnitude by *Doping*



Semiconductors

- SC resistivity between that of conductors and insulators
- Resistivity can be tailored over many orders of magnitude by *Doping*



Nanostructures

- Dimension(s) $< 100\text{nm}$

AlGaAs	$\sim 100\text{nm}$
GaAs	$\sim 20\text{nm}$
AlGaAs	$\sim 100\text{nm}$


Semiconductors

- SC resistivity between that of conductors and insulators
- Resistivity can be tailored over many orders of magnitude by *Doping*



Nanostructures

- Dimension(s) < 100nm



Spintronics




- Concerned with the generation, manipulation, and detection of spin polarization
- Technological example: HD read heads- “Spin Valve”
- Semiconductor Spintronics : No real world devices yet

How spin behaves in semiconductor material is currently being studied

Project Objectives

- Characterize electrical properties of semiconductor structures using the *Hall Effect*



$$n_s = IB / (qV_H)$$

n_s = Sheet Density

I = Current

B = Magnetic Field

q = Charge

V_H = Hall Voltage

$$\mu = 1 / (n_s e R_s)$$

R_s = Sheet Resistance

n_s = Sheet Density

e = Electron Charge

μ = Mobility


-
- Upgrade the PPMS (Physical Properties Measurement System) to allow van der Pauw measurements
 - Measure samples grown by MBE

System



- Current Source/DMM
- Control Computer
 - LabView for instrumentation control and data analysis
- PPMS Electronics
 - Temperature control
 - Magnetic field control
- Cryostat


Cryostat




Sample Puck



Resistivity



Hall



System



Current Source/DMM

Control Computer

LabView for instrumentation control
and data analysis


PPMS Electronics

Temperature control

Magnetic field control

Cryostat


Cryostat




Sample Puck



Resistivity



Hall



System



Current Source/DMM

Control Computer

LabView for instrumentation control
and data analysis


PPMS Electronics

Temperature control

Magnetic field control

Cryostat


Cryostat




Sample Puck



Resistivity



Hall



System



Current Source/DMM

Control Computer

LabView for instrumentation control
and data analysis


PPMS Electronics

Temperature control

Magnetic field control

Cryostat


Cryostat




Sample Puck



Resistivity



Hall



System



Current Source/DMM

Control Computer

LabView for instrumentation control
and data analysis


PPMS Electronics

Temperature control

Magnetic field control

Cryostat


Cryostat




Sample Puck




Resistivity



Hall




Sheet resistivity, sheet density, mobility




Spin lifetime depends strongly on carrier concentration

Sheet resistivity, sheet density, mobility




Spin lifetime depends strongly on carrier concentration

- Sheet Resistivity (Ohms/square)




Sheet resistivity, sheet density, mobility




Spin lifetime depends strongly on carrier concentration


- Sheet Resistivity (Ohms/square)



- Sheet Density (number/cm²)




Sheet resistivity, sheet density, mobility




Spin lifetime depends strongly on carrier concentration

- Sheet Resistivity (Ohms/square)



- Sheet Density (number/cm²)


- Mobility (cm²/V-s)



Calculating R_s , n_s , μ

I-V meas.


→ determine R → Calculate R_s



$$R_A = (R_{12} + R_{34})/2 \quad R_B = (R_{13} + R_{24})/2$$

$$\exp(-\pi R_A/R_S) + \exp(-\pi R_B/R_S) = 1$$


Calculate R_s



Calculating R_s , n_s , μ

I-V meas.


→ determine R → Calculate R_s




$$R_A = (R_{12} + R_{34})/2 \quad R_B = (R_{13} + R_{24})/2$$

$$\exp(-\pi R_A/R_S) + \exp(-\pi R_B/R_S) = 1$$

Calculate R_s




Hall measurement




$$n_s = (I/q) * (B/V_H)$$

$$= (I/q) * (1/\text{slope})$$


$$\mu = 1/(n_s e R_s)$$



Sheet resistivity, density, mobility vs. temperature




“multi-2DEG sample”



- Sheet Density

Generally increases with temp.




- Mobility

Strong function of impurities and temperature (phonons)

Remaining Tasks

- Sample puck modifications

Faster/Easier




- AC/Lockin measurement

Alternating current

More data

Signal/Noise



- Measure magnetic samples

Anomalous Hall Effect

“hysteresis”

